


ANTERIOR TEETH Special - Upper moulds

Mould type

- Aesthetic
- Youthful
- Prime



A2

40.3
9.5
7.3
(11)

A4

41.4
9.6
7.6
(7)

A7

43.5
9.9
7.9
(7)

A9

44.5
10.5
8.0
(8)

A11

48.0
11.1
8.8
(14)

A13

50.3
12.9
9.3
(16)

A43

43.2
9.2
7.9
(12)

A44

45.0
9.9
8.3
(13)

A52

45.8
10.4
8.3
(13)

A53

47.0
11.4
8.6
(16)

A64

47.0
9.7
8.6
(A71)

A72

49.0
10.4
9.1
(A73)

YM3

41.5
9.4
7.9
(11)

YM5

43.5
10.4
8.2
(9)

YW3

45.5
10.2
8.6
(13)

YW5

46.8
11.2
8.9
(14)

PM2

41.8
9.0
7.8
(11)

PM4

43.8
10.0
8.1
(12)

PM6

45.8
11.0
8.4
(13)

PM8

47.2
12.0
8.7
(14)

PW2

46.4
9.8
8.5
(13)

PW4

47.8
10.8
8.8
(A73)

PW6




48.5
11.8
9.1
(16)

PW8

50.0
12.8
9.2
(16)

ANTERIOR TEETH




Special - Lower moulds



POSTERIOR TEETH

Duratomic

Upper & Lower moulds



COMBINATION CHART

ANTERIORS		MOULD ARTICULATED WITH LOWER		COMBINES WITH POSTERIOR MOULD		
MOULD	A/6 mm	MOULD	A/6 mm	MOULD	UPPER P/4 mm	LOWER P/4 mm
A2	40.3	11	31.0	28	27.7	29.1
A4	41.4	7	32.0	28	27.7	29.1
A7	43.5	7	32.0	30	29.6	30.6
A9	44.5	8	33.1	32	31.6	33.0
A11	48.0	14	34.9	34	33.4	35.0
A13	50.3	16	36.7	34	33.4	35.0
A43	43.2	12	32.9	30	29.6	30.6
A44	45.0	13	34.7	32	31.6	33.0
A52	45.8	13	34.7	32	31.6	33.0
A53	47.0	16	36.7	34	33.4	35.0
A64	47.0	A71	35.0	32	31.6	33.0
A72	49.0	A73	36.5	34	33.4	35.0
YM3	41.5	11	31.0	30	29.6	30.6
YM5	43.5	8	33.1	32	31.6	33.0
YW3	45.5	13	34.7	32	31.6	33.0
YW5	46.8	14	34.9	32	31.6	33.0
PM2	41.8	11	31.0	30	29.6	30.6
PM4	43.8	12	32.9	30	29.6	30.6
PM6	45.8	13	34.7	32	31.6	33.0
PM8	47.2	14	34.9	34	33.4	35.0
PW2	46.4	13	34.7	32	31.6	33.0
PW4	47.8	A73	36.5	34	33.4	35.0
PW6	48.5	16	36.7	34	33.4	35.0
PW8	50.0	16	36.7	34	33.4	35.0

AVAILABLE SHADES

Vita® shades

A1 - A2 - A3 - A3.5 - A4 - B1 - B2 - B3 - B4 - C1 - C2 - C3 - C4 - D2 - D3 - D4

Myerson shade : B